

## EXPERIMENT 12M

FV 6/6/12

### INVESTIGATING FACTORS THAT AFFECT $\Delta T$ , $q_{\text{rxn}}$ , AND $\Delta H$ OF NEUTRALIZATION

**MATERIALS:** Coffee cup and lid, stir bar, magnetic stir plate, digital thermometer, two 100 mL graduated cylinders (for transport and measure of the acids and bases). Approximately 175 mL 1.0 M HCl, 50 mL 2.0 M HCl, 25 mL H<sub>2</sub>SO<sub>4</sub>, 300 mL 1.0 M NaOH, 50 mL 2.0M NaOH.

**PURPOSE:** To determine the temperature change, heat of reaction, and enthalpy change for a series of acid-base neutralization reactions.

**OBJECTIVES:** By the end of this experiment, the student should be able to:

1. Calculate  $\Delta H$  for a given reaction using thermodynamic data from tables.
2. Determine  $\Delta H$  experimentally for a given reaction using a coffee cup calorimeter.
3. Predict how thermodynamic values change when varying the quantity or concentration of reactants.

**PRE-LAB:** Read the entire lab guide and instructions and complete the attached pre-lab assignment.

### BACKGROUND:

The neutralization of an acid with a base in aqueous solution is an exothermic process. The heat transferred in this process ( $q_{\text{rxn}}$ ) can be measured using a calorimeter. Calorimeters can be very complicated and expensive (bomb calorimeters, isothermal calorimeters, differential scanning calorimeters, etc.) however a simple aqueous calorimeter can be made from any well insulated cup, such as a styrofoam coffee cup.

In this experiment, we will assume that no heat escapes or is absorbed by the coffee cup. Making this assumption, the energy transferred to the liquid inside the cup approximates the total energy released by the reaction. Therefore the heat of the calorimeter ( $q_{\text{cal}}$ ) approximates the heat transferred to the liquid ( $q_{\text{soln}}$ ) and it is the negative of  $q_{\text{rxn}}$ .

$$q_{\text{cal}} \approx q_{\text{soln}} \approx -q_{\text{rxn}} \quad (1)$$

The heat evolved by the neutralization reaction is transferred to the solution and can be measured as a temperature change where:

$$q_{\text{soln}} = m_{\text{soln}} \cdot s_{\text{soln}} \cdot \Delta T_{\text{soln}} \quad (2)$$

In addition, we will assume that the specific heat of the solution is the same as that of the solvent (water), therefore  $s_{\text{soln}} = 4.184 \text{ J/g } ^\circ\text{C}$ .

$\Delta H_{\text{rxn}}$  (the enthalpy change for a reaction) is often expressed in units of kJ/mol. When expressed in this unit,  $\Delta H_{\text{rxn}}$  is the amount of heat released or absorbed (at constant pressure) **per mole of reaction**. (The unit for  $\Delta H_{\text{rxn}}$  is sometimes written as kJ/mol<sub>rxn</sub> to emphasize this point.) In this lab you will calculate  $\Delta H_{\text{rxn}}$  as the enthalpy change associated with the formation of one mole of H<sub>2</sub>O through the process of acid/base neutralization. Since every mole of reaction produces one mole of water, this is also the enthalpy change for one mole of reaction. The total amount of heat produced by a reaction is given by:

$$q_{\text{rxn}} = n\Delta H_{\text{rxn}} \quad (3)$$

where n is moles of reaction (moles of water produced) so:

$$\Delta H_{\text{rxn}} = \frac{q_{\text{rxn}}}{n} = -\frac{q_{\text{cal}}}{n} \approx -\frac{q_{\text{soln}}}{n} \quad (4)$$

You will calculate  $q_{\text{soln}}$  using equation 2 and your measured temperature change. You will calculate the moles of water produced (n) using the moles of reactants added and the stoichiometry of the reaction. You will assume that the reaction goes to completion (the limiting reactant is completely consumed).

## PROCEDURE:

(**Caution:** Hydrochloric acid, sulfuric acid, and sodium hydroxide are corrosive and toxic.)

1. Construct a polystyrene-cup calorimeter, as demonstrated by your instructor.
2. For each of the five experimental conditions (below), perform the following protocol:
  - a. Place the acid in the calorimeter and monitor the temperature for 1 min, in order to ensure a stable initial temperature. Do not puncture the calorimeter with the thermometer.
  - b. Record the initial temperature.
  - c. Add the base to the calorimeter, and quickly stir the solution.  
(Make sure that you use a different graduated cylinder for the acid and base solutions to avoid mixing the reagents before the experiment is started.)
  - d. Monitor the temperature of the solution as it rises, and record the maximum temperature once it begins to drop again.
  - e. Dispose of the reaction solution by pouring down the drain and rinse the stir bar and calorimeter with distilled water. There is no need to dry the cup between expts.

**Experiment 1:** 50.0 mL 1.00 M hydrochloric acid and 50.0 mL 1.00 M sodium hydroxide.

**Experiment 2:** 100.0 mL 1.00 M hydrochloric acid and 100.0 mL 1.00 M sodium hydroxide.

**Experiment 3:** 50.0 mL 2.00 M hydrochloric acid and 50.0 mL 2.00 M sodium hydroxide.

**Experiment 4:** 25.0 mL 1.00 M hydrochloric acid and 75.0 mL 1.00 M sodium hydroxide.

**Experiment 5:** 25.0 mL 1.00 M sulfuric acid and 75.0 mL 1.00 M sodium hydroxide.

3. Calculate the  $\Delta T$  for each experiment. Compare your  $\Delta T$  values to those of the group next to you. If any  $\Delta T$  values differ by more than  $0.5^\circ\text{C}$ , repeat that experiment together.

## Clean-up/Disposal:

1. All solutions can be rinsed down the drain with plenty of water.
2. Rinse your calorimeter, stir bar, and any other glassware used.
3. Return equipment to their proper locations and clean up your lab area.

Name \_\_\_\_\_

Section # \_\_\_\_\_

**EXPERIMENTAL DATA – Expt. 12M:**

<u>Experiment #</u>	<u>Initial Temp (°C)</u>	<u>Maximum Temp (°C)</u>	<u>ΔT (°C)</u>
---------------------	--------------------------	--------------------------	----------------

**DATA TREATMENT:**

Calculate the following for each of your five experiments:

- number of moles of  $\text{H}^+$  initially present
- number of moles of  $\text{OH}^-$  initially present
- number of moles of  $\text{H}_2\text{O}$  produced
- temperature change
- heat released during the reaction (in joules)
- enthalpy for the reaction (in *kilojoules per mole of water formed*)

Record your answers in the table provided below:

<u>Experiment number</u>	<u>Initial mol <math>\text{H}^+</math>, mol</u>	<u>Initial mol <math>\text{OH}^-</math>, mol</u>	<u>mol <math>\text{H}_2\text{O}</math> produced, mol</u>	<u>ΔT, °C</u>	<u>Mass of solution g</u>	<u>q<sub>rxn</sub>, J</u>	<u>ΔH<sub>rxn</sub>, kJ/mol</u>
--------------------------	---	--	--	---------------	---------------------------	---------------------------	---------------------------------



- b. Experiments 1 and 4 each have a total reaction solution volume of 100.0 mL and use 1.00 M hydrochloric acid and 1.00 M sodium hydroxide. How do you account for their different  $q_{\text{rxn}}$  and  $\Delta T$ ?
- b. Experiments 4 and 5 each use 25.0 mL of a 1.00 M strong acid and 75.0 mL of 1.00 M sodium hydroxide. How do you account for their different  $q_{\text{rxn}}$  and  $\Delta T$ ?
4. Despite changes in reaction conditions (volume, concentration, monoprotic vs. diprotic acid),  $\Delta H_{\text{rxn}}$  (in kJ/mol) for each reaction was approximately the same. Explain why.
5. Based on a total reaction volume of 100.0 mL, what volumes of 1.00 M sulfuric acid and 1.00 M sodium hydroxide will produce the greatest amount of heat?



- d. Calculate the  $\Delta H$  (in kJ/mol) for the reaction forming 1 mole of water. (See equation 4.)
3. The neutralization of 50.0 mL of 1.00 M hydrochloric acid with 50.0 mL of 1.00 M sodium hydroxide causes a  $6.5^\circ\text{C}$  increase in temperature. Predict how the following changes to the experimental protocol would affect the value of the change in temperature. Explain your answers.

Hint: According to equations 1-3 given above:

$$\Delta T_{\text{soln}} = \frac{q_{\text{soln}}}{m_{\text{soln}} s_{\text{soln}}} = -\frac{q_{\text{rxn}}}{m_{\text{soln}} s_{\text{soln}}} = -\frac{n\Delta H_{\text{rxn}}}{m_{\text{soln}} s_{\text{soln}}}$$

Since  $\Delta H_{\text{rxn}}$  and  $s_{\text{soln}}$  are constants, the temperature change depends on the moles of water produced (moles of reaction) and the mass of the solution. Assume that the density of each solution is 1.0 g/mL.

- a. using 100.0 mL of 1.00 M hydrochloric acid and 100.0 mL of 1.00 M sodium hydroxide
- b. using 50.0 mL of 2.00 M hydrochloric acid and 50.0 mL of 2.00 M sodium hydroxide
- c. using 25.0 mL of 1.00 M hydrochloric acid and 75.0 mL of 1.00 M sodium hydroxide
- d. using 25.0 mL of 1.00 M sulfuric acid and 75.0 mL of 1.00 M sodium hydroxide (Be careful. Think about the value of n!)