

## Experiment 12C

2/27/2023

### THERMODYNAMIC PROPERTIES OF AMMONIUM CARBAMATE<sup>1</sup>

**MATERIALS:** Ammonium carbamate apparatus.

**PURPOSE:** The purpose of this experiment is to calculate the equilibrium constant,  $K_p$ , for a chemical reaction at several temperatures and to use these values to calculate the thermodynamic properties of the reaction.

**LEARNING OBJECTIVES:** By the end of this experiment, the student should be able to demonstrate the following proficiencies:


1. Derive an equation relating partial pressures of gaseous reactants and products to the equilibrium constant for a reaction.
2. Determine  $\Delta G^\circ$ ,  $\Delta H^\circ$  and  $\Delta S^\circ$  for a reaction from the values of the equilibrium constant at several temperatures.
3. Determine the  $\Delta G_f^\circ$  for a reactant from  $\Delta G^\circ$  for a reaction and the tabulated values of  $\Delta G_f^\circ$  for the products of the reaction.
4. Explain the correlation between  $\Delta G^\circ$  for a reaction and spontaneity under a given set of experimental conditions.

#### DISCUSSION:

The change in Gibbs free energy of a chemical reaction,  $\Delta G$ , is the criterion that is used to indicate whether or not the reaction is spontaneous.  $\Delta G$  is related to the enthalpy change of reaction,  $\Delta H$ , and the entropy change of reaction,  $\Delta S$ , by the expression  $\Delta G = \Delta H - T\Delta S$ , where  $T$  is the Kelvin temperature. Knowledge of the thermodynamic quantities associated with a chemical reaction is important since experimental conditions, especially temperature, can sometimes be manipulated to cause the reaction to occur in the desired direction.

For example, in a closed atmosphere, such as in submarines and space vehicles, control of the carbon dioxide level is crucial. Normal air contains a very small quantity of carbon dioxide, approximately 0.03% by volume. But high levels of carbon dioxide are toxic. The average person produces 120 g of  $\text{CO}_2$  per day at normal activity levels. This amount increases as the activity level increases. High carbon dioxide concentrations, measured as the partial pressure of  $\text{CO}_2$  ( $P_{\text{CO}_2}$ ), can cause impaired judgment, dizziness, loss of muscular coordination, unconsciousness, and eventually death. For this reason, methods of removing carbon dioxide are of great interest to the military. For more information on carbon dioxide toxicity and removal methods, consult the *Naval Applications of Chemistry* website, at the link for submarine air treatment.<sup>2</sup>

Early NASA spacecraft used non-regenerative carbon dioxide scrubbers such as lithium hydroxide. They remove carbon dioxide by the spontaneous chemical reaction indicated by the following equation:




Sufficient lithium hydroxide must be carried on board the spacecraft to absorb the carbon dioxide produced during the entire flight. Assuming 100% efficiency of scrubbing, 130 g or 5.5 moles of lithium hydroxide are needed per day per crew member.

---

<sup>1</sup> Joncich, M.J., Solka, B.H., Bower, J.E., *J. Chem. Ed.*, **1967**, *44*, 598. Bennett, R.N., et al., *Trans. Faraday Soc.*, **1953**, *49*, 925.


<sup>2</sup> [http://intranet.usna.edu/ChemDept/\\_files/documents/navapps/ADDITIONAL-INFO/Submarine%20Air/SUBMARINE%20AIR%20TREATMENT\\_2010.pdf](http://intranet.usna.edu/ChemDept/_files/documents/navapps/ADDITIONAL-INFO/Submarine%20Air/SUBMARINE%20AIR%20TREATMENT_2010.pdf)

On a submarine, which must stay submerged for extended periods of time, carrying enough lithium hydroxide to remove all of the carbon dioxide produced by the crew would require too much space and weigh too much<sup>3</sup>. Instead, submarines use a regenerative scrubbing system. This regenerative system is based upon the chemical reaction between monoethanolamine (MEA) and carbon dioxide shown in the equation below.




The position of equilibrium for this reaction can be controlled by changing the temperature. At low temperatures and relatively high carbon dioxide pressures, carbon dioxide is absorbed because the forward reaction is spontaneous. At higher temperatures and low carbon dioxide pressures, carbon dioxide is produced because the reverse reaction is spontaneous. An exchange system is used which allows for removal of between 70% and 90% of the carbon dioxide from the crew air spaces. The waste carbon dioxide is collected and pumped overboard.


The purpose of this experiment is to study a system similar to the MEA/CO<sub>2</sub> system, the ammonium carbamate system. Ammonia, like MEA, reacts with carbon dioxide to form a carbamic acid (H<sub>2</sub>NCO<sub>2</sub>H). The carbamic acid then reacts with a second ammonia molecule to form a carbamate salt, ammonium carbamate (NH<sub>4</sub>CO<sub>2</sub>NH<sub>2</sub>).



The overall reaction is:



In this experiment the thermodynamic properties of the reverse reaction, the decomposition of ammonium carbamate, will be determined by measuring the equilibrium vapor pressure as a function of temperature. Ammonium carbamate decomposes to form carbon dioxide and ammonia as shown in the equation below.



The thermodynamic properties of a reaction can be determined from equilibrium constants for the reaction at several temperatures. These values are related by the equation

$$\ln K_p = - \frac{\Delta H^\circ}{R} \frac{1}{T} + \frac{\Delta S^\circ}{R} \quad (7)$$

where  $K_p$  is the equilibrium constant. This assumes that the values of  $\Delta H^\circ$  and  $\Delta S^\circ$  remain essentially constant over the entire temperature range. Since this equation has the form of the equation for a straight line,  $y = mx + b$ , the slope and intercept of a plot of  $\ln K_p$  versus  $1/T$  can be used to calculate values for  $\Delta H^\circ$  and  $\Delta S^\circ$ . The desired thermodynamic information for a reaction and its components can be calculated from this equation, other basic equations of thermodynamics, and experimental data.

## PROCEDURE:

See pre-lab. Experimental data must be collected **the week before** the lab is scheduled. A diagram and picture of the apparatus are shown on the next page.

---

<sup>3</sup>For a three month cruise, the submarine would have to carry about 13 kilograms of lithium hydroxide per person or 1.4 tons for a 100 person crew. This would require about 602 L or 23 ft<sup>3</sup> of space.

Name \_\_\_\_\_


Section \_\_\_\_\_

Partner \_\_\_\_\_

Date \_\_\_\_\_

**DATA SECTION**  
**Experiment 12C**

<b>Table 1 - Data for the Reaction: <math>\text{NH}_4\text{CO}_2\text{NH}_2 (\text{s}) \rightleftharpoons 2 \text{NH}_3 (\text{g}) + \text{CO}_2 (\text{g})</math></b>			
Reading	Temperature (°C)	Hg height (cm) high	Hg height (cm) low
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			



**Figure 1. Figure and picture of apparatus.** The closed-end mercury manometer containing ammonium carbamate is completely immersed in a temperature-controlled bath. There is a meter stick for measuring the heights of the mercury, and a thermometer for measuring the temperature of the bath (and manometer). Hg “high” and Hg “low” refer to the measured position of the mercury meniscus on the “left” (reference,  $P = 0 \text{ atm}$ ) and right” (sample) sides of the manometer, respectively.



### Part B: Calculating Thermodynamic Functions

1. Using the results from your  $\ln K_P$  vs.  $1/T$  graph, calculate the value for  $\Delta H^\circ$ . Show your work. It is assumed that over the small temperature range of this experiment,  $\Delta H^\circ$  is independent of temperature. Calculation of  $\Delta H^\circ$  from the slope of your line is therefore valid.  $R = 8.314 \text{ J/mol}\cdot\text{K}$

$$\Delta H^\circ = \text{_____ kJ/mol}$$

2. Using the values of the slope and the y-intercept, calculate the value of  $K_P$  at 298.15 K.

$$K_{P, 298.15} = \text{_____}$$

3. Calculate the value of  $\Delta G^\circ_{298.15}$  using your value for the equilibrium constant  $K_P$  at 298.15 K.

$$\Delta G^\circ_{298.15} = \text{_____ kJ/mol}$$

4. Using your experimentally determined values of  $\Delta H^\circ$  and  $\Delta G^\circ_{298.15}$ , calculate  $\Delta S^\circ$ . Note that the value of  $\Delta S^\circ$  is also assumed to be constant over the experimental temperature range.

$$\Delta S^\circ = \text{_____ J/K}\cdot\text{mol}$$

5. Calculate the change in the internal energy for the reaction,  $\Delta E^\circ$ , at 298.15 K from your value of  $\Delta H^\circ$ .

$$\Delta E^\circ = \Delta_r H^\circ - (\Delta n)RT \quad \text{where } \Delta n = \text{the change in the total moles of gas (from balanced chemical equation).}$$

$$\Delta E^\circ = \text{_____ kJ/mol}$$

6. The thermodynamic quantities calculated to this point ( $\Delta H^\circ$ ,  $\Delta G^\circ_{298.15}$ ,  $\Delta S^\circ$  and  $\Delta E^\circ$ ) are measures of the changes in the enthalpy, Gibbs free energy, entropy, and internal energy, respectively, at 298.15 K and standard state for the overall reaction. Values of the standard enthalpies of formation ( $\Delta H_f^\circ$ ), standard free energies of formation ( $\Delta G_f^\circ$ ), and the standard molar entropies ( $S^\circ$ ) for  $\text{NH}_3(\text{g})$  and  $\text{CO}_2(\text{g})$  are likewise tabulated in various sources, including your textbook and in the *Handbook of Chemistry and Physics*. By combining the overall thermodynamic reaction quantities with the specific quantities for  $\text{NH}_3$  and  $\text{CO}_2$ , we can calculate values for the standard heat of formation ( $\Delta H_f^\circ$ ), standard free energy of formation ( $\Delta G_f^\circ$ ), and the standard molar entropy ( $S^\circ$ ) for ammonium carbamate.
- a. Thermodynamic quantities for  $\text{NH}_3(\text{g})$  and  $\text{CO}_2(\text{g})$  at 298.15 K are tabulated in Appendix 4, Table A4.3 on pp. APP-18 thru APP-24 of the Gilbert text. Fill in the table below, using all significant figures from these sources.

**Data Treatment Table 1**

Compound	$\Delta H_f^\circ_{298.15}$ (kJ/mol)	$\Delta G_f^\circ_{298.15}$ (kJ/mol)	$S^\circ_{298.15}$ (J/mol·K)
$\text{NH}_3(\text{g})$			
$\text{CO}_2(\text{g})$			

Use the thermodynamics values you calculated for the *overall reaction* ( $\Delta H^\circ$ ,  $\Delta G^\circ_{298.15}$ ,  $\Delta S^\circ$ ) and the values from Data Treatment Table 1 (values for the *individual* components of  $\text{NH}_3$  and  $\text{CO}_2$ ) to determine the thermodynamic quantities for *ammonium carbamate alone* in questions b, c, and d below.

- b. Calculate the value of  $\Delta H_f^\circ$  for  $\text{NH}_4\text{CO}_2\text{NH}_2$  at 298.15 K.

$$\Delta H_f^\circ = \text{_____ kJ/mol}$$

c. Calculate the value of  $\Delta G_f^\circ$  for  $\text{NH}_4\text{CO}_2\text{NH}_2$  at 298.15 K.

$\Delta G_f^\circ =$  \_\_\_\_\_ kJ/mol

d. Calculate the value of  $S^\circ$  for  $\text{NH}_4\text{CO}_2\text{NH}_2$  at 298.15 K.

$S^\circ =$  \_\_\_\_\_ J/K $\oplus$ mol





Name \_\_\_\_\_

Section \_\_\_\_\_

Date \_\_\_\_\_

**PRE-LAB EXERCISES**  
**Experiment 12C**


Complete these exercises by the due date specified by your instructor.

**Part A: Data Collection**

1. You must collect the data for this experiment **the week before** the lab is scheduled. The apparatus required for this experiment will be set up daily (Monday – Friday) on the chemistry lab deck between the hours of 0745 and 1530. The temperature of the apparatus will be changed in the morning (before 0745) and at lunch (after 1200) so it is possible to obtain 2 measurements each day.
2. Using the apparatus, carefully **record the temperature** and **the heights of the mercury levels in the manometer** in Table 1 of the DATA SECTION, using the correct number of significant figures and units. It is your responsibility to complete Table 1 by the date set by your instructor. Alternatively, your instructor may provide a Google spreadsheet to collect the class data.

**Part B: Derivation of Equations**

1. Write the equilibrium expression,  $K_P$ , for the decomposition reaction studied in this experiment.



2. In the experimental apparatus, all the gases inside the sealed end were formed from the decomposition of solid ammonium carbamate, as shown above. Applying Dalton's Law of Partial Pressure, what is the expression for the total pressure in the sealed end in terms of the gases present?

$$P_{\text{total}} = \text{_____} + \text{_____}$$

3. a. Based on the balanced chemical equation, what is the stoichiometric mole ratio between  $\text{NH}_3$  and  $\text{CO}_2$ ?

$$\text{NH}_3 : \text{CO}_2$$

$$\text{_____} : \text{_____}$$

- b. Using Dalton's Law of Partial Pressures and the stoichiometric relationship between ammonia and carbon dioxide (and therefore between  $P_{\text{NH}_3}$  and  $P_{\text{CO}_2}$ ), what fraction of the total pressure in the tube comes from carbon dioxide? What fraction of the total pressure in the tube comes from ammonia?

$$\text{NH}_3: \text{_____} \quad \text{CO}_2: \text{_____}$$

4. Substitute the relationships just derived above into the equation for  $K_P$  so that  $K_P$  is given in terms of only the total pressure,  $P_{\text{total}}$ , rather than the partial pressures of ammonia and carbon dioxide. This equation is needed on page 4, #3.