

## Experiment 4E

FV 6/20/16

### THE DECOMPOSITION OF POTASSIUM CHLORATE

**MATERIALS:** Two test tubes: (18x150), clamp, ring stand, Bunsen burner, weighing boat, glass wool, pure  $\text{KClO}_3$ , unknown mixture containing  $\text{KClO}_3$  and inert substance,  $\text{MnO}_2$


**PURPOSE:** The purpose of this experiment is to study the decomposition of potassium chlorate, both by verifying the identity of one of the products and by quantitatively determining the correct stoichiometry.

**LEARNING OBJECTIVES:** By the end of this experiment, the student should be able to demonstrate the following proficiencies:

1. Explain the relationship between the mass of a substance and the number of moles of a substance.
2. Apply stoichiometric ratios between the moles of reactant(s) and product(s) in balanced chemical reactions.
3. Calculate the percent by mass of a compound in a mixture.
4. Explain the purpose of a catalyst.

#### DISCUSSION:

*Stoichiometry.* A major emphasis of chemistry is the understanding chemical reactions. This requires knowing the correct formulas for all reactants and products involved in the reaction, as well as the relative molar amounts of each. Such information is provided by the balanced chemical reaction, but where does *that* come from? The answer is that reactions are determined by experiment. Careful mass measurements and physical and/or chemical tests allow one to deduce the proper reaction. Only when that is understood can one start to consider useful applications of the reaction. Consider the title reaction, the thermal decomposition of potassium chlorate. When  $\text{KClO}_3$  is heated strongly, it breaks down releasing oxygen gas and leaving behind a thermally stable (i.e., heat-insensitive) solid residue of an ionic potassium compound.



There are at least three plausible reactions one can write for the process, but only one occurs to any significant extent. Which one is actually observed can only be determined by experiment, such as those conducted here. By measuring the amount of oxygen lost when a sample of potassium chlorate is heated, we will be able to determine the stoichiometric coefficients of  $\text{KClO}_3$  and  $\text{O}_2$  in the reaction, and thus determine the correct reaction.


*Relevant Naval Application.* On submarines, oxygen for breathing is normally produced through electrolysis of water. Details relating to this process will be studied later in the course. In an emergency, a chemical process is used to produce oxygen gas for breathing, specifically the decomposition of sodium chlorate at high temperature (i.e., above  $300^\circ\text{C}$ ), producing oxygen gas and a solid sodium salt. Unfortunately, there are several complications associated with this reaction which must be remedied if the production of oxygen gas for breathing is to be performed safely and efficiently in this practical application.

First, though the decomposition reaction occurs at temperatures above  $300^\circ\text{C}$ , it is extremely slow and therefore impractical for oxygen production in bulk. This is remedied by adding a *catalyst*, in this case manganese (IV) oxide, which significantly increases the rate of the reaction, without itself being consumed.

Second, the intense flame used to raise the temperature of the sodium chlorate above  $300^\circ\text{C}$  is produced by a combustion reaction, which *consumes* large quantities of oxygen gas, whereas the purpose of the overall process is to *produce* oxygen gas. While this issue cannot be completely remedied, small amounts of iron metal are mixed in, reacting with some of the oxygen to produce iron oxide and releasing large quantities of energy which helps maintain the mixture above the  $300^\circ\text{C}$  decomposition temperature. After the candle is ignited, the oxygen-consuming flame used to initiate the decomposition reaction is replaced by this iron combustion process, making it more self-sustaining.

Third, while the desired decomposition reaction predominates, there is another decomposition reaction which produces toxic chlorine gas, oxygen gas and sodium oxide. This is remedied by including small amounts of barium peroxide in the mixture, which reacts with the toxic chlorine gas to produce barium chloride and oxygen gas.

In summary, the “chlorate” or “oxygen” candle used for emergency production of oxygen gas for breathing on submarines consists of a mixture of sodium chlorate, iron, a small amount of barium peroxide, and a fibrous binding material. In practice, each candle burns near 400°C for 45-60 minutes, and produces approximately 115 SCF (standard cubic feet) of oxygen gas at 0.5 psig (pounds per square inch, gauge pressure), which is enough oxygen for about 100 people. The stored candles represent a significant fire hazard since they are self-sustaining in oxygen.



*Use of potassium chlorate.* In this experiment, potassium chlorate will be used instead of the sodium chlorate employed commercially (see Figure 1). As you should suspect, analogous reactions occur, with all of the same complications. The only remedy that will be applied here will be the inclusion of the manganese (IV) oxide catalyst. Since all of the procedures will be carried out in the fume hood, any toxic chlorine gas produced will be safely carried away in the ventilation system. Why is NaClO<sub>3</sub> used commercially, rather than KClO<sub>3</sub>? The principal reason is cost; sodium salts are typically much less expensive than their potassium counterparts.

*Energetics.* As should be clear from the discussion of the Navy’s chlorate candles, practical applications almost always have to consider energy changes associated with reactions. Most chemical and physical processes are accompanied by changes in energy – some release energy as they proceed, and some require an input of energy in order to sustain the process. We will also examine some elementary concepts of energy changes associated with the reactions observed.

*Material Safety Data Sheets and International Chemical Safety Cards.* Any institution where chemicals are used is required to have copies of the material safety data sheets (MSDS) available for use. These sheets provide key information relating to health hazards, appropriate storage, handling and disposal arrangements, fire and explosive hazards, required control measures, physical/chemical properties, and reactivity data. In this experiment, the MSDS for potassium chlorate will be used to help guide the experimental study of its decomposition reactions. In general, prior to any chemical procedure, the relevant MSDS should be consulted to assure safe and proper procedures are followed. Another system which provides similar information is the International Chemical Safety Card system. Both MSDS and Safety Cards are available on-line through links found on the Plebe Chemistry homepage.

## PROCEDURE:

### Part A: Determination of the Stoichiometry of the Decomposition Reaction of Potassium Chlorate

1. Add a small amount of manganese (IV) oxide (about half as large as a pea) to a clean, dry test tube and determine the mass of the tube and manganese (IV) oxide on an *analytical* balance. Record this mass, using the proper number of significant figures and units, in the Data Section.
2. Use the *top-loading balance* to pre-weigh about 1.5 g of pure potassium chlorate into a plastic weighing boat. It is not necessary to record this mass yet. If any clumps of potassium chloride are present, break them up with a spatula, and then add the potassium chlorate to the test tube containing the MnO<sub>2</sub>. **Be careful not to introduce any foreign material into the bottles of potassium chlorate since explosive mixtures could be produced.**
3. Determine the mass of the test tube and its contents on the same analytical balance used initially (in step A.1) and record the mass in the Data Section.
4. Carefully tap the test tube against your hand to mix the contents to obtain a uniformly gray mixture. Do NOT insert any stirring devices, to avoid removing material and changing the mass! Tapping the test tube should also knock down to the bottom any materials which may be adhering to the upper inside walls.
5. Place a *loose* plug of glass wool in the mouth of the test tube. This will allow oxygen to escape but prevent any solid from spilling out of the tube. Determine the total mass of the tube and its contents (including glass wool) on the same analytical balance used previously and record the mass in the Data Section.
6. Clamp the test tube to a ring stand at a slight angle as shown in Figure 2. Be sure that the clamp does not have plastic sleeves as these will burn during the experiment. Place the clamp near the open end of the test tube so that the clamp will not melt while the test tube is being heated. (Don't squeeze the clamp too tightly as you may crack the tube.) Be sure the open end of the test tube is not pointed toward anyone or toward the lab aisles.
7. **Carefully observe what happens as you heat the sample.** Heat the tube gently at first since oxygen is driven off quickly as the decomposition of the potassium chlorate begins. Move the burner around to achieve uniform heating. Increase the rate of heating as the rate of gas evolution decreases, finally heating as strongly as possible for an additional three or four minutes.
8. Allow the tube to cool to room temperature and determine the mass of the tube and its total contents, including the residue and the glass wool, on the same analytical balance used previously. Record the mass in the Data Section.
9. If time permits, heat the same sample a second time and record its mass when cooled.

☞ Answer Questions #1 - #3 on page E4E-9.




Figure 2. Experimental set-up.

**Part B: Determination of the Percent by Mass of Potassium Chlorate in an Unknown Mixture**

1. From your instructor, obtain an unknown sample containing a mixture of potassium chlorate and an inert substance. **Record the unknown number in the space provided in the Data Section.**
2. Repeat the procedure outlined in Part A, using approximately 1.5 g of your assigned unknown in place of the pure potassium chlorate. Record all masses, using the proper number of significant figures and units, in the Data Section.

**Clean Up:**

1. Place all used test tubes, including their contents, in the **designated solid waste container** in the laboratory, NOT in the broken glass box.

Name \_\_\_\_\_

Section \_\_\_\_\_

Partner \_\_\_\_\_

Date \_\_\_\_\_

**DATA SECTION**  
**Experiment 4E**

**Part A: Determination of the Stoichiometry of the Decomposition Reaction of Potassium Chlorate**

Mass of tube + MnO<sub>2</sub> \_\_\_\_\_

Mass of tube + MnO<sub>2</sub> + KClO<sub>3</sub> \_\_\_\_\_

Mass of KClO<sub>3</sub> \_\_\_\_\_

Mass of tube and total contents (including glass wool)

    Before heating \_\_\_\_\_

    After heating \_\_\_\_\_ \*

Mass of oxygen gas evolved \_\_\_\_\_ \*

**Part B: Determination of the Percent by Mass of Potassium Chlorate in an Unknown Mixture**

Unknown number \_\_\_\_\_

Mass of tube + MnO<sub>2</sub> \_\_\_\_\_

Mass of tube + MnO<sub>2</sub> + unknown \_\_\_\_\_

Mass of unknown \_\_\_\_\_

Mass of tube and total contents (including glass wool)

    Before heating \_\_\_\_\_

    After heating \_\_\_\_\_ \*

Mass of oxygen gas evolved \_\_\_\_\_ \*

\* If time permits, heat the sample a second time.

**DATA TREATMENT**  
**Experiment 4E**

**Part A. Determination of the Stoichiometry of the Decomposition Reaction of Potassium Chlorate**

In this section the balanced chemical equation for the decomposition of potassium chlorate will be determined. From the mass information recorded in the Data Section, you will calculate the number of moles of potassium chlorate present in the test tube before heating (question A.1) and the moles of oxygen gas generated (question A.2). After these two calculations, determine the stoichiometric ratio between potassium chlorate and oxygen gas (A.3) and generate the correct balanced chemical equation (A.4). This equation should also include the chemical formula for the solid material left in the test tube after heating. *Show all work.*


(A.1) From the mass of potassium chlorate used, calculate the number of moles of potassium chlorate consumed.

(A.2) From your data, calculate the number of moles of oxygen gas (O<sub>2</sub>) evolved.

(A.3) Based on your answers above, calculate the numerical value of the following mole ratio (your instructor may ask you to record your ratio on the board):

$$\frac{\text{moles } O_2 \text{ produced}}{\text{moles } KClO_3 \text{ decomposed}} =$$

(A.4) Balance the three reactions below, as done in the Pre-Lab, and show the ratio moles O<sub>2</sub>/moles KClO<sub>3</sub> for each reaction. Which of the three reactions best fits your experimental data? Explain your answer below.



### Part B. Determination of the Percent by Mass of Potassium Chlorate in an Unknown Mixture

Having determined the balanced chemical equation for the decomposition of potassium chlorate (from page E4E-6), it is now possible to determine the percent by mass of potassium chlorate in an unknown mixture containing  $\text{KClO}_3$  and an inert ingredient. As before, heating the sample will cause the decomposition to occur and oxygen gas will be produced. By calculating the number of moles of oxygen produced (B.1), the moles and mass of potassium chlorate originally present in the unknown sample (B.2) can be determined. From that the mass percent potassium chlorate in the unknown can be obtained (B.3). Show your work.

(B.1) From your data, calculate the number of moles of  $\text{O}_2$  evolved when the unknown was heated.

(B.1) Using the correct stoichiometric equation, calculate the mass of potassium chlorate decomposed to provide the  $\text{O}_2$  loss observed.

(B.1) From your data, calculate the percent by mass of potassium chlorate in the unknown sample. Pay attention to significant figures.

$$\text{mass \%} = \frac{\text{mass KClO}_3}{\text{mass unknown}} \times 100$$

Unknown number: \_\_\_\_\_

Percent by mass of  $\text{KClO}_3$  in unknown: \_\_\_\_\_

**QUESTIONS FOR CONSIDERATION: (Your instructor will determine which questions you must answer.):**

1. Ignoring any side reactions and assuming the reaction occurs completely, how large (in kg) an oxygen candle ( $\text{KClO}_3$ ) would be needed to supply 8 people with enough oxygen for 24 hours on a small submarine? Although this depends on the size of the person and their respiration rate (activity), according to NASA<sup>1</sup>, an average person needs about 0.84 kg of  $\text{O}_2$  per day.
2. Starting with 100.0 g of potassium chlorate, and assuming that 95% decomposes into potassium chloride/oxygen and 5% decomposes into potassium oxide/chlorine/oxygen, how many grams of residue would remain after complete reaction? If the investigator had assumed that only the potassium chloride/oxygen reaction was taking place, what would he find for a mole ratio of oxygen gas: potassium chlorate?
3. In practice, rather than using a Bunsen burner to maintain the heat needed for this reaction, a small amount of iron metal is introduced into the mixture. Write the balanced reaction with oxygen, if iron (II) oxide is the product. Explain why this reaction is only useful if it produces large quantities of energy.
4. If the test tube had not been fully cooled off before weighing, buoyance effects in the balance tend to make the mass appear lighter than it actually is. In this case, what would be the effect on the calculated mole ratio of oxygen : potassium chlorate?
5. When used for producing oxygen, the small percentage of decomposition which follows the second pathway produces hazardous chemicals. In practice, small amounts of barium peroxide are introduced to react with these products, as described in the Discussion section above. Balance this reaction.
6. What would have been the expected mass of residue if the potassium chlorate reaction produced potassium hypochlorite instead of potassium chloride, assuming no alternate decomposition pathways? Determine the oxygen to potassium chlorate mole ratio for this scenario. Does this possible reaction pathway appear to contribute at all, based on the experimental data?
7. If experimental errors cause the number of moles of oxygen determined in the data analysis to be about 3% low, will an incorrect chemical equation for the decomposition reaction be obtained? Explain your answer; some numbers will help prove your point.

---

<sup>1</sup>Wieland, P.O., Designing for Human Presence in Space: An Introduction to Environment Control and Life Support Systems, NASA Reference Publication 1324, 1994, pp. 6, 183-262.



Name \_\_\_\_\_

Section \_\_\_\_\_

Date \_\_\_\_\_



**IN-LAB QUESTIONS**  
**Experiment 4E**



Complete these questions during lab.

1.
  - a. Record your observations as the mixture is heated.
  
  
  
  
  
  
  
  
  
  
  - b. How do you know when the decomposition of potassium chlorate is near completion?
  
  
  
  
  
  
  
  
  
  
2.
  - a. What is the purpose of  $\text{MnO}_2$  in this experiment?
  
  
  
  
  
  
  
  
  
  
  - b. Why doesn't the  $\text{MnO}_2$  affect the stoichiometry of the reaction?
  
  
  
  
  
  
  
  
  
  
3. Why do you think it would be advisable to heat the sample a second time?

Name \_\_\_\_\_

Section \_\_\_\_\_

**PRE-LAB QUESTIONS**  
**Experiment 4E**

1. Select which statements(s) are true, according to the Material Safety Data Sheet for potassium chlorate.  
(<https://www.fishersci.com/shop/msdsproxy?productName=AAA170750I&productDescription=potassium-chlorate-&catNo=AAA170750I&vendorId=VN00024248&storeId=10652>)

- a)  $\text{KClO}_3$  is a liquid at room temperature and 1 atm.
- b) Ingestion may cause hemolysis, methemoglobinemia, cyanosis, anuria, coma, and convulsions.
- c)  $\text{KClO}_3$  is a powerful oxidizing agent
- d) Contact with combustible substances may cause fire or explosion

2. Since this experiment involves high temperatures, what are the melting point and decomposition temperature (both in  $^{\circ}\text{C}$ ) for potassium chlorate?

Melting Temperature: \_\_\_\_\_

Decomposition Temperature: \_\_\_\_\_

\*NOTE: Some MSDSs list the decomposition temperature as the “boiling” temperature.

3. Based on the values in problem 2, what will you see happening to the potassium chlorate solid as you heat it to high temperatures?

- a) the potassium chlorate will decompose before it melts
- b) the potassium chlorate will melt before it begins to decompose
- c) the potassium chlorate will melt and decompose simultaneously

4. Following are three possible reactions that could occur when a sample of  $\text{KClO}_3$  (s) is decomposed. One goal of this lab is to determine the correct one. Balance all three reactions. (Use the lowest possible integer coefficients.)

- a) \_\_\_  $\text{KClO}_3$ (s)  $\rightarrow$  \_\_\_  $\text{KClO}_2$ (s) + \_\_\_  $\text{O}_2$ (g)
- b) \_\_\_  $\text{KClO}_3$ (s)  $\rightarrow$  \_\_\_  $\text{KClO}$ (s) + \_\_\_  $\text{O}_2$ (g)
- c) \_\_\_  $\text{KClO}_3$ (s)  $\rightarrow$  \_\_\_  $\text{KCl}$ (s) + \_\_\_  $\text{O}_2$ (g)

You might wish to copy your results to question A.4 on p. E4E-6 as you will need them during lab.

5. Using the correct stoichiometric coefficients from the balanced equations in question 4, write the ratio: moles  $\text{O}_2$  produced/ moles  $\text{KClO}_3$ (s) decomposed for each reaction:

Reaction 4(a)\_\_\_\_\_

Reaction 4(b)\_\_\_\_\_

Reaction 4(c) \_\_\_\_\_

6.  $\text{NaClO}_3$  is used on Navy submarines for what purpose? (NOTE: The chemistry of  $\text{NaClO}_3$  is very similar to that of  $\text{KClO}_3$ . However,  $\text{NaClO}_3$  has the advantage in that it is less expensive and  $\text{KClO}_3$ ).

- a)  $\text{NaClO}_3$  candles are an emergency source of light
- b)  $\text{NaClO}_3$  is used as a source of oxygen for emergency situations
- c)  $\text{NaClO}_3$  is added to paint to prevent corrosion
- d)  $\text{NaClO}_3$  is used as a food additive.